

AVANT™ Stealth Audiometer

Features
Integrated
Amplifiers
For Free Field
Testing

Clinical Audiometry Has Reached A New Level Of Performance

This dual channel audiometer allows fast and accurate pure tone, speech and bone testing with an intuitive user interface for data presentation and patient monitoring.

Separate outputs for all transducers allow them to remain plugged in and calibrated adding speed and convenience to the testing protocol.

The Stealth is PC-based and includes your choice of standard transducers (TDH 39 or Insert Earphones). This system has powerful internal amplifiers and is NOAH™, Sycle.Net & TIMS® Compatible.

Options include high frequency and Quick SIN™ testing capability.

MedRx®



AD-S-ADH138-1

www.medrx-usa.com or Call Toll Free: 888-392-1234