

# HANDS FREE!

## Motorized Polar HIT<sup>+</sup> Chamber



Small Footprint  
w: 11" d: 8.5" h: 5.2"

The New AVANT™ Polar HIT<sup>+</sup> Chamber Conducts A Battery Of Automated, Selectable ANSI & IEC Tests.

This HIT Box Also Provides A Motor Driven Rotating Test Pad That Automatically Generates Directional Polar Plots For Fast, Accurate Testing.

# MedRx®

[www.medrx-usa.com](http://www.medrx-usa.com) or Call Toll Free: 888-392-1234

# AVANT Polar HIT<sup>+</sup>

## Motorized Hearing Instrument Test Chamber



### AVANT Polar HIT<sup>+</sup> Performs The Following Measurements:

- Directionality
- OSPL 90
- Full On Gain
- Reference Test Gain
- Frequency Response
- Equivalent Input Noise
- Harmonic Distortion
- Battery Drain
- Input - Output
- Attack / Release
- Induction Coil



OT-S-ADH131-1

### Additional Features:

- Omni and Directional Microphone Verification
- New Automatic Rotating Test Pad Provides Convenient Directional Hearing Instrument Testing
- Measures Polar Plots Of Instrument Sensitivity
- Allows For Customized Test Protocols
- Simple One-Cable USB Connectivity Eliminates The Need For External Power Supply
- Superior Acoustic Isolation For Accurate Testing
- Lighted Visual Indicators Show When A Test Is Active
- Results Are Stored In NOAH™ And Easily Printed

**MedRx**INC®

[www.medrx-usa.com](http://www.medrx-usa.com)

1200 Starkey Rd. #105,

Largo, FL 33771

Toll Free 888-392-1234

727-584-9600 • Fax: 727-584-9602

MedRx's Authorized Representative in Europe  
(Regulatory affairs only) Emergo Europe • Molenstraat 15  
2513 BH The Hague • The Netherlands  
Tel: (31) (0) 70 345-8570 • Fax: (31) (0) 70 346-7299