

# AVANT™ Polar HIT<sup>+</sup>

**Automated  
Polar Plot  
Directional  
Microphone  
Testing**



## Motorized Hearing Instrument Test Chamber

This unique system is designed to test virtually all hearing instruments with the additional automated capability of testing directional instruments. The standard battery of ANSI and IEC tests are selectable and conducted automatically by the system. The test chamber connects to a USB port and is NOAH™, Sycle.Net & TIMS® Compatible.

### AVANT™ Polar HIT<sup>+</sup> Performs The Following Measurements:

- Directionality • OSPL 90 • Full On Gain
- Reference Test Gain • Frequency Response
- Equivalent Input Noise • Input - Output
- Harmonic Distortion • Battery Drain
- Attack / Release • Induction Coil

**MedRx**®

OT-S-ADH139-1



[www.medrx-usa.com](http://www.medrx-usa.com) or Call Toll Free: 888-392-1234