

# The **UltraVac**

## Hearing Instrument Restoration System



Featuring  
Pulse  
Power

The UltraVac™ hearing instrument restoration system, has many features including Pulse Power. The vacuum wand will pull up to 28 inches of mercury when the constant suction feature is utilized. The Pulse Power feature greatly increases suction by pulsing the vacuum in the wand. Should wax become lodged in the vacuum needle, the user may simply switch on the pressure feature and reverse the airflow. Once the switch is engaged, the drying chamber vacuum will run for 4 minutes and shut off automatically.

*MedRx also offers an optional  
hard-side carrying case*

UV-S-ADH138-1

# MedRx®

**www.medrx-usa.com or Call Toll Free: 888-392-1234**