

# Battery Operated Light Source

## Install Freshly Charged Batteries Every Day

A Second Set Of 3 Batteries Is Supplied  
So That A Fully Re-Charged Set Is Always Available




Battery Charger



### Charging

Main Features:

- Micro Computer Control
- AA, AAA Battery Types are compatible
- Suitable for Ni-CD and Ni-MH rechargeable batteries
- Four charge channels that can charge one to four batteries
- Folding AC plug blade for easy storage



Type	Size	CAP.	Charging Time	
			1,2 cells	1,2,3,4 cells
NI-MH	AA	2000mA	Approx. 120min	Approx. 240min
	AAA	600mA	Approx. 80min	Approx. 150min
NI-CD	AA	2500mA	Approx. 150min	Approx. 300min
	AAA	850mA	Approx. 120min	Approx. 200min



Insert Battery Pack into the Light Source Handle with "UP" at Top of Battery Pack

